

User Manual

Version 3.0

Product name: I-5™ 2X High-Fidelity Master Mix

Cat #: I5HM-100, I5HM-200, I5HM-5000, I5HM-25k, I5HM-OEM

Description:

I-5TM High-Fidelity DNA Polymerase is an ultra-high fidelity and high processivity enzyme. It produces the most accurate copies of DNA, and performs at an ultra-high rate. I-5TM is perfect for applications like cloning, in vitro amplifing materials for protein expression, SNP analysis by sequencing, and high-specificity PCR.

Convenient and Usage:

I-5TM 2X High-Fidelity Master Mix has a buffer system that is designed to offer maximum performance of the I-5TM High-Fidelity DNA Polymerase. It offers the highest fidelity and incredibly robust performance. The convenient master mix allows reaction to be set up in room temperature and only requires the additions of primers and DNA template for the amplification.

Application:

- High-specificity PCR amplification
- High-throughput PCR
- Various cloning technologies
- Difficult amplification

Recommended Storage Condition: -20°C

Instructions:

	25 μl Reaction	50 μl Reaction	Final Concentration
I-5™ 2X High-Fidelity Master Mix	12.5 µl	25 μΙ	1X (see notes)
10 μM Primer A	1 µl	2 μΙ	400 μM
10 μM Primer B	1 µl	2 μΙ	400 μM
Template DNA	as required	as required	see notes
Water (ddH ₂ O)	up to 25 μl	up to 50 µl	

Thermocycling Conditions

3 Step PCR

Step	Temperature	Time	
Initial Denaturation	98°C	2 minutes	
Denaturation	98°C	10 seconds	25-35 cycles
Annealing	45°C - 68°C	10-15 seconds	
Extension	72°C	15-30 seconds / kb	
Final Extension	72°C	1-5 minutes	
	4°C	Hold	

2 Step PCR (see notes)

Step	Temperature	Time	
Initial Denaturation	98°C	2 minutes	
Denaturation	98°C	10-15 seconds	25-35 cycles
Combined Annealing & Extension	72°C	15-30 seconds / kb	
Final Extension	72°C	1-5 minutes	
	4°C	Hold	

Notes:

Recommended DNA Template addition

Genomic DNA 50-250 ng
Plasmid DNA 1 pg-10 ng
Viral DNA 1 pg-10 ng

2 Step PCR

Use of 2 Step PCR is recommended when the primer's Tm values are $>68^{\circ}$ C

Mg⁴

The final concentration of Mg^{2+} in 1X I-5TM PCR Master Mix is 2 mM. Add additional Mg^{2+} as needed in 0.5 mM increments

Suggested the final Mg²⁺ concentration from 2 mM to 4 mM

PCR Product / Cloning

The 2X I- 5^{TM} PCR Master Mix results in PCR product with blunt ends Blunt-end cloning is recommended after PCR with this product

For T/A-cloning, the PCR product should be purified before A-addition