

# User Manual

Version 3.0

**Product name:** T4 DNA Ligase

**Cat #:** TL-100, TL-200, TL-300, TL-400, TL-OEM, B-TD10

## Description:

T4 DNA Ligase catalyzes the formation of a phosphodiester bond between the terminal 5'-phosphate and 3'-hydroxyl groups of duplex DNA or RNA. The enzyme efficiently joins blunt and cohesive ends and repairs single-stranded nicks in duplex DNA, RNA or DNA/RNA hybrids.

## Protocol:

Set up the following reaction in a microcentrifuge tube on ice. The T4 DNA Ligase Buffer should be thawed and resuspended at room temperature.

(T4 DNA Ligase should be added last. Note that the table shows a ligation using a molar ratio of 1:3 vector to insert.)

Components: (20  $\mu$ l reaction)  
10X T4 DNA Ligase Buffer\* 2  $\mu$ l  
Vector DNA (3 kb) 50 ng (0.025 pmol)  
Insert DNA (1 kb) 50 ng (0.076 pmol)  
Nuclease-free water to 20  $\mu$ l  
T4 DNA Ligase 1  $\mu$ l\*

- Gently mix the reaction by pipetting up and down and microfuge briefly.
- For cohesive (sticky) ends, incubate at 16°C overnight or room temperature for 10 minutes.
- For blunt ends or single base overhangs, incubate at 16°C overnight or room temperature for 2 hours (alternatively, high concentration of T4 DNA Ligase can be used in a 10 minute ligation).
- Chill on ice and transform 1-5  $\mu$ l of the reaction into 50  $\mu$ l competent cells.