



Molecular Cloning Laboratories

Cell Line Identification

Sample Shipping Instructions

Carefully follow the instructions below to insure your cell line identification samples are successfully received by our laboratory. For questions, please contact one of our specialists at 1-650-872-0245.

Step 1

If sending cells, prepare a cell pellet containing at least 10,000 cells.

Step 2

Ship the cell pellet on blue ice packs or dry ice.

Step 3

If performing your own DNA isolation, please send 10 - 20 μ l of purified DNA (at a minimum concentration of 1 ng/ μ l). We recommend using TE-4 buffer (10mM Tris-HCl, 0.1mM EDTA [pH 8.0]). Please provide the concentration of DNA if known.

Step 4

DNA samples in solution should be shipped on blue ice packs or dry ice. Dried DNA samples can be shipped at room temperature.

Step 5

Fill out the shipment form and send the samples to our laboratory via overnight courier. Shipping address:

MCLAB (Attn: Cell Line)
320 Harbor Way
South San Francisco, CA 94080
USA