

User Manual

Version 2.0
Revision Date: 09/21/2015

Product name: YT Agar Plates

Cat #: YTA-100, YTA-200, YT-AC501, YT-AC502, 2YT-Amp101, 2YT-Amp102, 2YT-ACa501, 2YT-ACa502, 2YT-Ach341, 2YT-Ach342, 2YTA-100, 2YTA-200, 2YT-AK301, 2YT-AK302

Description:

0.8% Tryptone, 0.5% yeast extract, 0.5% NaCl, 0.7% agar

Recommended storage condition: 4°C

Name	Description	Application
YT Agar Plates	0.8% Tryptone, 0.5% yeast extract, 0.5% NaCl, 0.7% agar	For use with <i>E. coli</i> and other enteric bacteria.
2YT Agar Plates	1.6% tryptone, 1.0% yeast extract, 0.5% NaCl, 1.5% agar, 100µg/ml ampicillin, 20µg/ml tetracycline, 2.0% glucose	For growth and propagation of <i>E.coli</i> infected with the single strand filamentous bacteriophage M13.