



Material Safety Data Sheet

Version 3.0

Revision Date: 06/03/2011

Print Date: 06/03/2011

Product and Company Identification

Product Name: Taq DNA Polymerase (regular)
Product Use: For Research Use Only. Not for use in diagnostic procedures
Company: Molecular Cloning Laboratories(MCLAB)
320 Harbor Way
South San Francisco
CA 94080
Telephone: 888-MCLAB-88
650-871-8771
Fax: 650-871-8796

Composition/Information on Ingredients

Hazardous/Non-hazardous Compounds

CHEMICAL NAME	CAS #	Formula	SARA 313
Glycerol	56-81-5	HOCH ₂ CHOHCH ₂ OH	NO

Synonyms: glycerin, glycerol USP, glycerine, 1,2,3-propanetriol, propanetriol, 1,2,3-trihydroxypropane, bulbold, citifluor AF 2, cristal, emergy 916, glyrol, glycerol ophthalgan, glyciterol, glycy alcohol, osmoglyn, pricerine 9091

Hazardous Ingredients Identification

Generally not hazardous in normal handling, however good laboratory practices should always be used. Avoid long term exposure to skin or by inhalation.

First Aid Measures

FIRST AID:

SKIN: Wash exposed area with soap and water. If irritation persists, seek medical attention.

EYES: Wash eyes with plenty of water for at least 15 minutes, lifting lids occasionally. Seek Medical Aid.

INHALATION: Remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen.

INGESTION: Give several glasses of milk or water. Vomiting may occur spontaneously, but it is not necessary to induce. Never give anything by mouth to an unconscious person.

Fire Fighting Measures

Fire Extinguisher Type:

Carbon Dioxide, dry chemical powder or appropriate foam

Fire/Explosion Hazards:

Slight Flammability. Material must be heated before ignition will occur.

Fire Fighting Procedure:

Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and clothing.

Accidental Release Measures

Absorb spill with inert material, then place in a chemical waste container. Dispose of in a manner consistent with federal, local law.

Handling and Storage

HANDLING	STORAGE	NOTES:
Avoid inhalation	Keep tightly closed	Hygroscopic
Avoid contact with eyes, skin, and clothing	Store in cool, dry place	
Avoid prolonged or repeated exposure		

Exposure Controls/Personal Protective Equipment

Engineering:

Safety Shower and eye wash. Ensure adequate ventilation.

Personal Protective Equipment:

Respiratory:

NIOSH/MSHA-approved respirator

Hand:

Compatible chemical resistant gloves

Eye:

Compatible safety goggles

Skin and Body:

Lightweight protective clothing

Hygiene Measures:

Accordance with good industrial hygiene and safety practices

Environmental Exposure Controls:

Prevent product from entering drains

Exposure Limits:

CHEMICAL NAME	OSHA PEL (TWA)	OSHA PEL (Ceiling)	ACGIH OEL (TWA)	ACGIH OEL (STEL)
Glycerol	15 mg/m ³ total dust 5 mg/m ³ respirable fraction	-	10 mg/m ³	-

Physical and Chemical Properties

General:

Physical Form: Liquid

Appearance: Clear

Odor: Faint

Physical Properties:

Melting Point: 18°C

Boiling Point: 289°C

Flash Point: 193°C

Autoignition Temp: 400°C

Explosion Limits in Air:

Lower: 0.9%

Specific Gravity: 1.2610g/cm³

Solubility: Miscible in water

Vapor Pressure: 0.0025 mm Hg @ 5°C

Vapor Density: 3.17 g/L

pH: 5.5-8.0

Stability and Reactivity

Stability:

Stable

Materials to Avoid:

Strong oxidizing agents, strong bases

Hazardous Decomposition products:

Carbon Monoxide, Carbon Dioxide

Hazardous Polymerization:

Will not occur

Toxicological Information

Conditions aggravated / target organs: Single dose oral toxicity is low. Ingestion of larger amounts may cause injury. Signs and symptoms of excessive exposure may be central nervous system effects and increased blood sugar levels. Person with preexisting conditions may be more susceptible.

DOT Classification: Not Regulated

DOT regulations may change from time to time. Please consult the most recent version of the relevant regulations.

Ecological Information

Data not available.

Disposal considerations

Dispose of in accordance with local regulations.

Transport Information

Substance considered non-hazardous for ground or air transport. No DOT shipping name required.

Regulatory Information

International Inventories:

CHEMICAL NAME	TSCA	PICCS	ENCS	DSL	NDSL	AICS
Glycerol	Listed	Listed	Listed	Listed	-	Listed

U.S. Federal Regulations:

SARA 313: Not regulated

U.S . State Regulations:

CHEMICAL NAME	MA - RTK
Glycerol	Listed

WHMIS Hazard Class:

Product has been classified in accordance with the hazard criteria of the CPR and the MSDS contains all of the information required by the CPR.

Other Information

Disclaimer: This product intended for R&D and clinical diagnostic uses. Use of material must be supervised by technically qualified individual. This above information is believed to be correct but may not be all inclusive and shall be used as a reference guide only. Company shall not be held liable for any damage resulting from handling or contact with above product.

End of MSDS