

# User Manual

Version 1.0  
Revision Date: 07/30/2013

**Product name:** DEPC-Treated Water

**Cat #:** DTM-100, DTM-200, DTM-300, DTM-400

## Description:

DEPC-treated Water is deionized, high quality, molecular biology grade water.

## Recommended Storage Conditions:

Room temperature.

## Protocol:

DEPC-Treated Water is suitable for use with RNA without any DNase or RNase activity detected. It is prepared by incubating with 0.1% diethylpyrocarbonate (DEPC), autoclaved and sterile filtered.